

## Relion 3.1 $\beta$ -galactosidase tutorial



#### Last Updated: 2021-03-10

Leicester Institute of Structural and Chemical Biology



#### Processing CryoEM Data

Recent innovations in camera technology and processing capabilities have resulted in the so-called 'resolution revolution' in cryo-electron microscopy. While we are now able to calculate cryoEM maps at very high resolution, the field is still young. As a result, the level of experience available in most structural biology labs for single-particle cryoEM analysis (SPA) is not yet on-par with that of NMR or Xray crystallography. This tutorial is meant to provide some of that knowledge in a structured way.

This tutorial uses Relion (<u>J Struct Biol, 2012 Dec 180(3). 519-530</u>, <u>github.com/3dem/relion</u>) for almost all processing tasks, and is based very heavily on the Relion tutorial (<u>ftp.mrc-Imb.cam.ac.uk/pub/scheres/relion31\_tutorial.pdf</u>) written by Sjors Scheres and co-workers. If you've never done SPA before – complete the entire tutorial before moving to other packages. While the details of how each software package works will change, the underlying concepts won't.

In addition to the Relion tutorial, we have collected much of our SPA experience here. As our experience and expertise evolves, we will hope to capture more information in this document. For example, we now routinely use AI based autopicking via Relion's External job type. While this is not covered in the tutorial, there is a section at the end describing how to do this. We have also included specific parameters for how to submit jobs to our scheduling system, including differences between the tutorial data and an actual project.

It's important to emphasize that this is a collaborative production. We have incorporated suggestions and critiques from PI's, facility managers, Post-docs, Ph.D. students and undergraduate students. If you feel you have something to contribute to this document, please let us know at <u>cryoem@leicester.ac.uk</u>.

Finally: about the style of this tutorial. SPA is done almost entirely on computer clusters that use the Linux operating system. If you don't know Linux, or how to use the command-line, you can still use this tutorial – but you may find you capabilities greatly expanded by a separate Linux tutorial. Some of the parameters we use are extremely precise. If you have a problem, please check the parameters you're using very carefully before raising an issue. The tutorial dataset is small, but you should be able to produce a reasonable quality map with it. If you are not able to successfully produce a map better than 4 Å – go back, and try again. You should not attempt to process real data until you can successfully process this data.



# Downloading files from the image processing cluster (Windows)

Some parts of this tutorial require you to download maps to your computer. Here's how to do it if you have a Windows machine.

- 1. Install WinSCP (from *Program Installer* or <u>winscp.net</u>)
- 2. Open WinSCP
- 3. Add a new Site (Session > New Session...)
  - 1. Make sure protocol is 'SFTP'
  - 2. Fill in
    - Host name: 143.210.183.163
    - User name: your cluster username
    - · Cluster password
  - 3. Setup Tunnel
    - 1. Click Advanced
    - 2. Choose the Tunnel page
    - 3. Fill in:
      - Host name: spectre2.le.ac.uk
      - User name
      - · University password

New Site	Session File protocol: SFTP V		Environment     Directories     Recycle bin     Encryption     SFTP     Shell	Connect through SSH tunnel Host to setup tunnel on Host name: spectre2.le.ac.uk	Port number:
	Host name: 143.210.183.163 User name: Pr YYYY Save V	Port number: 22 • assword: Advanced	Connection Proxy SSH Key exchange Authentication Bugs Note	User name: Passwo [x0x] Tunnel options Local tunnel port: Tunnel authentication parameters Private key file:	Autoselect v
Tools 🔻	Manage 🔻 🔁 Login 🔻	Close Help	Color V	ОК	Cancel Help



# Downloading files from the image processing cluster (Mac)

Some parts of this tutorial require you to download maps to your computer. Here's how to do it if you have a Mac.

- 1. You will need a 'tunnel' to get the data:
  - Open a terminal window (Applications > Utilities > Terminal), and run: ssh -D 2222 -q -C -N XXX@spectre2.le.ac.uk
     Replacing XXX with your UoL username. You'll need to enter your UoL password. If the command succeeds, it not generate any errors and appear to hang. This is normal.
- 2. Download FileZilla Client (filezilla-project.org)
- 3. Uncompress the download by double-clicking on the .bz2 file
- 4. Open FileZilla
- 5. Setup Generic Proxy
  - 1. Open Settings (Edit > Settings...)
  - 2. Select Generic Proxy page
  - 3. Fill in
    - Type: Socks 5
    - Host: localhost
    - Port: 2222
  - 4. Click 'OK'

	Settings
Select page:	Generic proxy
<ul> <li>Connection</li> <li>FTP</li> <li>Active mode</li> <li>Passive mode</li> <li>FTP Proxy</li> <li>SFTP</li> <li>Generic prov</li> <li>Transfers</li> <li>FTP: File Types</li> <li>File exists action</li> <li>Interface</li> <li>Passwords</li> <li>Themes</li> <li>Date/time format</li> </ul>	Type of generic proxy: None HTTP/1.1 using CONNECT method SOCKS 4 SOCKS 5 Proxy host: Proxy port: Proxy password: Note: Using a generic proxy forces passive mode on FTP connections.

- 6. Add a new 'Site'
  - 1. Open Site Manager (File > Site Manager...)
  - 2. Click on 'New site'
  - 3. Change the name
  - 4. Fill in
    - Protocol: SFTP
    - Host:143.210.183.163
    - User: your Username
  - 5. Click 'Connect'
  - 6. Enter your cluster password
  - 7. Agree to 'Trust Host?'

	agei
Select entry:	General Advanced Transfer Settings Charset
My Sites	Protocol: SFTP - SSH File Transfer Protocol
	Host: 143.210.183.163
	Logon Type: Ask for password
	User: tjr22
	Password:
	Background color: None
$\frown$	Comments:
2 New site New folder	
New Bookmark Rename	
Delete Duplicate	
	Connect OK Cancel



#### Accessing the image processing cluster with NoMachine

NoMachine is the easiest way to access the image processing cluster. While not as powerful as other methods, it's more than enough for this tutorial. NoMachine works with Windows, Mac, and Linux.

- University Windows computer Install using *Program Installer* 
  - NoMachine
  - UCSF Chimera
- Non-Univsersity Windows machine, Mac, Linux
  - 1. Install NoMachine from: www.nomachine.com
  - 2. Install UCSF Chimera from: www.cgl.ucsf.edu/chimera
  - 3. Download the UoL HPC connection file for NoMachine from:

uniofleicester.sharepoint.com/sites/Research-Computing/SitePages/nomachine-spectre.aspx

4. Install the UoL connection into NoMachine by double-clicking on it.



#### Accessing the image processing cluster with NoMachine

NoMachine is the easiest way to access the image processing cluster. While not as powerful as other methods, it's more than enough for this tutorial. NoMachine works with Windows, Mac, and Linux.

- Logging in to the cluster
  - 1. Double-click on the UoL HPC connection, then enter your UoL username and password Hover mouse right here
  - 2. Once logged in, change your view settings by hovering the mouse over the upper right-hand corner of the window (you may see the window slide across as you move to the corner)
  - Click on the 'Display' icon 3.
  - Click on the 'Change settings' icon 4.
  - 5. Set your Resolution to at least 900 high (e.g. 1440x900) Resolution 800x480
  - 6. Click Done, then on the 'Fit to window' icon, then Done
  - 7. Open a terminal, and log-in to the image processing cluster with: ssh –Y XXX@143.210.183.163

Where XXX is your UoL username. You'll be prompted for your password. Use your cluster password, not your UoL password. If you are asked if you want to trust this host, answer yes

#### Click here to open Terminal



## 🎹 🚰 Tue 23 Mar, 15:3 Click to view month calendar 2 service OK







#### Setting up tutorial project

Once logged in (make sure to connect X11), setup the Relion directory structure using the following commands:

1. Make the Relion project top-directory and a directory to contain the micrograph movies (commands are case-sensitive)

```
mkdir relion_tutorial
cd relion_tutorial
mkdir Micrographs
```

- 2. Symlink (DON'T COPY) the micrograph movies, and the gain reference, into your Micrographs directory (include the trailing . in your commands.)
  - cd Micrographs
    ln -s /net/common/relion\_3.1\_tutorial/Micrographs/\*.tiff .
    ln -s /net/common/relion\_3.1\_tutorial/Micrographs/gain.mrc .
    cd ..
- 3. Make sure that you've symlinked the data by listing the directory (you should see several tiff files and the gain reference file listed.)

ls Micrographs

3. Examine the acquisition notes

```
cat /net/common/relion_3.1_tutorial/Micrographs/NOTES
```

4. Start Relion (and return your prompt).

relion &

• Relion will ask you if you want to start a new project directory here. Answer 'Yes'.





#### Importing micrographs / micrograph movies

Relion keeps track of what we do, in order to simplify the 'pipeline' of data processing for us. The first step is to import our data.

	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial	To start a new Relion Job, click on the job type in
ile Jobs Schedules	Movies/mics Others Running	
Import Motion correction	Import raw movies/micrographs? Yes 🗘 ?	the Available Jobs box.
CTF estimation Manual nicking	Raw input files: Micrographs/*.tiff ?	
Auto-picking	Are these multi-frame movies? Yes 🗘 ?	
Particle extraction		1. CIICK ON Import
2D classification	Optics group name: opticsGroup1 ?	
3D initial model	MTF of the detector: mon/MTF/mtf_k2_300kV.star ? Browse	
3D classification 3D auto-refine	Pixel size (Angstrom): 0.885	2. On Movies/mics tab:
3D multi-body	Voltage (kV): 200 200 200 200 200 200 200 200 200 20	
Bavesian polishing	Spherical aberration (mm): 1.4 7	
Mask creation	Amplitude contrast: 0.1	Import raw movies/micrographs?
Join star files Particle subtraction	Beamtilt in X (mrad): 0	
Post-processing	Beamtilt in Y (mrad): 0	Yes
External		
Literitor	Schedule Check command Continue!	
		Raw Input files:
I/O view Job actions	Current: 001: Import/job001/ Display:	Migrographs / t tiff
Finished jobs	Running jobs Input to this job	MICIOGIaphs/ CIII
	001: Import/iob001/	
		Are these Multi-frame movies.
		Are chese Marci-Irane movies.
		Yes
	Scheduled jobs Output from this job	
		MTF of the detector:
		/net/common/MTF/mtf k2 300kV.star
importing		
Written Import/job001/movie done!	es.star with 0 items (0 new items)	
		Pixel Size (Angstrom):
		(From NOWER)
		(FIOIII NOTES)
		Voltago (kV).
		Voicage (KV).
		(From NOTES)
		Spherical aberration (mm):
		(From NOTES)
		(1.611 1.6126)
		Amplitude contrast:
		U • 1
		3. UICK RUN! UNCE

You can examine information Relion provides information about many of the options by clicking on the small ?. Examine the information for Amplitude contrast before proceeding.

This job should run very quickly. In the output window you can see how many items were imported (should be 24.) Note that there's a small idiosyncrasy that sometimes prevents Relion's GUI from updating a job's status unless you click on either it or on another job. Click on the import job now and it should move to the Finished jobs box.



#### **STAR Files**

Relion stores metadata in .star files. You can look at these files to find out more details about your processing. Let's examine the import star file in the terminal

```
ls
cd Import
ls
cd job001
ls
less movies.star
          # version 30001
          data optics
          loop
          _rlnOpticsGroupName #1
          rlnOpticsGroup #2
          _rlnMtfFileName #3
          rlnMicrographOriginalPixelSize #4
          _rlnVoltage #5
          rlnSphericalAberration #6
          rlnAmplitudeContrast #7
          opticsGroup1
                        1 ../../net/common/MTF/mtf_k2_300kV.star
                                                                              0.885000
          200.000000 1.400000
                                    0.100000
          # version 30001
          data_movies
          loop_
          rlnMicrographMovieName #1
          _rlnOpticsGroup #2
          Micrographs/20170629_00021_frameImage.tiff
                                                             1
          Micrographs/20170629 00022 frameImage.tiff
                                                             1
```

#### This is a basic STAR file. It has the following structure

- One data block (data\_optics), called 'optics'
  - One table in that data block (loop\_), without a name
  - Seven columns in the table, the first of which is called 'rInOpticsGroupName' (\_rlnOpticsGroupName)
  - One row in the table.
- A second data block called 'movies'
  - etc,...

Press 'q' to quit less. Type cd ~/relion\_tutorial to return to the relion\_tutorial directory within your home directory.



#### **Motion Correction (Relion's implementation)**

In order to reduce the impact of beam-induced motion on our images, we record movies. The first step in processing is to align the frames of those movies and then sum them up into a single image. Here, we'll use Relion's implementation of the MotionCor2 algorithm. We don't normally use dose-weighting.

	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Motion Running
Import	Input movies STAR file: Import/job001/movies.star   Browse
CTF estimation	First frame for corrected sum: 1
Auto-picking	Last frame for corrected sum: -1
Particle extraction Subset selection	Dose per frame (e/A2): 1.277 7
2D classification	Pre-exposure (e/A2): 0
3D initial model 3D classification	Do doso weighting? No
3D auto-refine	Save pop-dose weighted as well? No
3D multi-body CTF refinement	Save sum of power spectra? No
Bayesian polishing	Sum power spectra every e/A2: 4
Mask creation Join star files	
Particle subtraction	
Local resolution	
External	
	Schedule Check command KUN!
I/O view Job actions	Current: Give_alias_here Display:

	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Motion Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D startification	Bfactor: 150 7 Number of patches X, Y 1 1 7 Group frames: 1 7 Binning factor: 1 7 Gain-reference image: Micrographs/gain.mrc 7 8 Browse
2D classification 3D initial model 3D classification 3D auto-refine 3D multi-body CTF refinement	Gain rotation: No rotation (0)
Bayesian polishing Mask creation Join star files Particle subtraction Post-processing Local resolution	Use RELION's own implementation? No
I/O view Job actions	Schedule         Check command         Run!           Current:         Give_alias_here         Display:         +

X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial File Jobs Schedules I/O Motion Running Import Motion correct Number of MPI procs: 1 Idelian correction Manual picking Auto-picking Auto-picking Particle extraction 2D classification 3D initial model 3D classification 3D cl Number of threads: 1 Submit to queue? Yes \$ ? Queue name: gpu Queue submit command: sbatch 1 Maximum run time (d-h:m:s) 3-00:00:00 GPU count (0-4) Standard submission script: //net/prog/slurm/relion.sh Minimum dedicated cores per node: ? Browse - 7 Additional arguments: 2 External Schedule I/O view Job actions Current: Give\_alias\_here Display:

- 1. Click on the Motion Correction job in the Available jobs box
- 2. On I/O tab:

Input movies Star file: Import/job001/movies.star

Dose per frame: (From NOTES)

Do dose-weighting? <mark>No</mark>

3. On Motion tab:

Number of patches X,Y 5 5

Gain-reference image: Micrographs/gain.mrc

Use RELION's own implementation? Yes

4. On Running tab:

Number of MPI processes: 1 Number of Threads:

4 (should be divisor of number of frames)

Submit to Queue?: <mark>Yes</mark>

GPU count:

Click Run! ONCE

This job should take  $4\sim5$  minutes. On real datasets, use 4 MPI, but make sure that MPI × threads  $\leq 32$  (this is true for all jobs.) To ensure the best performance, make sure you run the job exclusively by changing sbatch to sbatch --exclusive.



#### **CTF Estimation with CtfFind4**

Once we have motion-corrected micrographs, we can estimate the CTF parameters which will be required later on for reconstruction. We prefer CtfFind4 (via the Relion GUI) for this.

	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O CTFFIND-4.1 Gctf Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification 3D indication 3D auto-refine 3D auto-refine Bayesian polsing May Star files Particle subtraction Post-processing Local resolution	Input micrographs STAR file: 2/corrected_micrographs.star ? Browse Use micrograph without dose-weighting? Ino ? ? Estimate phase shifts? No ? ? Phase shift - Min, Max, Step (deg) 0 180 10 ? Amount of astigmatism (A): 100 ?
External	Schedule Check command Run!
I/O view Job actions	Current: Give_alias_here Display:

00	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O CTFFIND-4.1 Gctf Running
Import Motion correction CFF estimation Manual picking Particle extraction Subset selection 2D classification 3D initial model 3D classification 3D classification CFF refinement Bayesian polishing Mask creation Join star files Particle subtraction Post-processing Local resolution	Use CTFFIND-4.1? Yes ¢ ? CTFFIND-4.1 executable: /net/prog/ctfind4/ctffind Use power spectra from MotionCorr job? No ¢ ? Use exhaustive search? No ¢ ? Estimate CTF on window size (pix) 1 7 FFT box size (pix) 512 7 Minimum resolution (A): 50 7 Maximum defocus value (A): 5000 7 Maximum defocus value (A): 5000 7 Defocus step size (A): 500 7 Defocus step size (A): 500 7 C
External	Schedule Check command Run!
I/O view Job actions	Current: Give_alias_here Display:

1. On I/O tab:

Input micrographs STAR file: MotionCorr/job002/corrected\_micrograph s.star

Use micrograph without dose weighting: No (Only say <u>yes</u> if you did dose-weighting and also saved non-dose-weighted images)

Estimate phase shifts? No (Only say Yes for phase-plate data)

2. On CTFFIND-4.1 tab: Use CTFFIND4.1? Yes FFT box size (pix): 512 (Use 1024 or 2048 for datasets that have very weak signals)

	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
ile Jobs Schedules	I/O CTFFIND-4.1 Gctf Running
Import Motion correction	Number of MPI procs: 3
Manual nicking	Submit to queue? Yes \$ ?
Auto-picking	Queue name: gpu ?
Particle extraction	Queue submit command: sbatch
Subset selection	Maximum run time (d-h:m:s) 3-00:00:00 ?
3D initial model	GPU count (0-4) 0
3D classification	Standard submission script: /net/prog/slurm/relion sh
3D auto-refine 3D multi-body	Minimum dedicated cores per node: 1
CTF refinement	
Bayesian polishing	Additional arguments:
Ioin star files	
Particle subtraction	
Post-processing	
Local resolution	
External	Schedule Check command Run!

3. On Running tab:

Number of MPI procs: 3 (20 for real data) Submit to queue? Yes GPU count (0-4): 0 Click Run! ONCE

This job should take 10~30 seconds. On real datasets, use 16 MPI.



#### Selecting good micrographs via CTF estimation

The first step in finding good particles is to select good micrographs. The easiest way to do this is by eliminating micrographs with obviously bad CTF estimates. We'll use Relion's subset selection job for this task

			XR	ELION-3.1-beta-c	ommit-7f8	d1d: /home/tji	22/relion31-t	tutorial		
File	Jobs	Schedules	I/O	Class options	Subsets	Duplicates	Running			
Impo Motio CTF e Manu Auto Partie 2D c 3D in 3D c 3D a 3D m CTF r Baye Mask Join s Partie Post- Loca	ort con corre estimat ual pick pickin cicle extr et sele lassifica uto-refin uto-refin uto-refin uto-refinem sian po creati star file cle sub proces l resolu	action ing action action ation odel ation ne dy dy traction sing traction sing tion		Select cla OR select fr OR sele OR selec	sses from rom microg cct from pa ct from picl	model.star: graphs.star: (ticles.star: (ked coords:	l/job003/mic	rographs_ctf.star	200	Browse Browse Browse Browse
I/O	view	Job actions	Cu	rrent: Give ali	as here	(	Schedule Disp	Check comman	a (	Run!
<u> </u>										



1. On I/O tab:

OR select from micrographs.star: CtfFind/job003/micrographs\_ctf.star (Make sure all the other fields are empty.)

2. Alias the job In Current: all\_good\_micrographs

Click Run!

On Relion Display GUI:

- 1. Click the box beside Sort images on:
- 2. Sort on rlnDefocusU
- 3. Nr. Columns: 3
- Click Display!

Each image has three parts. The background image is the power spectrum from your micrograph. The upper-left hand quadrant is the circularly averaged power spectrum, after taking into account the fitted astigmatism. The wedge in the upper right-hand quadrant is the fitted power spectrum.



- 1. Right-click and choose 'Invert selection', this will select all micrographs.
- Left click on any micrographs that don't have a good match between the calculated and measured CTF. Pay close attention to the first few and last few micrographs
- 3. Right-click and select Save STAR with selected images.
- Repeat the selection procedure, sorting on rlnDefocusU, rlnDefocusV, rlnCtfFigureOfMerit,

rlnCtfAstigmatism, and

- rlnCtfMaxResolution.
  - 1. Close the CTF display window.
  - 2. Change the Sort images on: selection
  - 3. Display!



### Manual picking

The best way to assess the quality of your data is via manual picking. We can also use the picked particles to generate 2D class averages for reference-based autopicking. Note here that the Particle diameter refers only to the size of the circle Relion draws on the micrographs and has nothing to do with the extracted particle box size used.

	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Display Colors Running
Import Motion correction CTF estimation Manual picking Particle extraction Subset selection 20 classification 30 initial model 30 classification 30 auto-critica	Input micrographs: <u>plect/job004/micrographs.star</u> <b>8 Browse</b>

	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Display Colors Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification 3D initial model 3D auto-refine 3D auto-refine 3D auto-refine Bayesian polishing Mask creation	Particle diameter (A): 20 + 1 + 1 + 20 + 1 + 20 + 1 + 20 + 1 + 20 + 1 + 20 + 1 + 20 + 1 + 20 + 1 + 20 + 1 + 20 + 20
Particle subtraction Post-processing	Scale for CTF image: 1

1. On I/O tab:

Input micrographs: Ctffind/job003/micrographs\_ctf.star

2. On Display tab:

Particle diameter (A): 20

Click **Run!** ONCE. This will open the Manualpicking GUI

	Contraction Contra						
Fil	e						
	aphs/20170629_00021_frameImage.mrc	pick	0	CTF	10869.1	Ĥ	
	aphs/20170629_00022_frameImage.mrc	pick	0	CTF	9752.89		
	aphs/20170629_00023_frameImage.mrc	pick	0	CTF	10673.5		
	aphs/20170629_00024_frameImage.mrc	pick	0	CTF	11676.6		
	aphs/20170629_00025_frameImage.mrc	pick	0	CTF	10673.5		
	aphs/20170629_00026_frameImage.mrc	pick	0	CTF	8351.08		
	aphs/20170629_00027_frameImage.mrc	pick	0	CTF	9148.28		
	aphs/20170629_00028_frameImage.mrc	pick	0	CTF	10126.8	) U	
	aphs/20170629_00029_frameImage.mrc	pick	0	CTF	8276.56		
	aphs/20170629_00030_frameImage.mrc	pick	0	CTF	9056.51		
	aphs/20170629_00031_frameImage.mrc	pick	0	CTF	10115.9		
	aphs/20170629_00035_frameImage.mrc	pick	0	CTF	9304.42		
	aphs/20170629_00036_frameImage.mrc	pick	0	CTF	10056.7		
	aphs/20170629_00037_frameImage_mrc	nick		CTE	11063.9		

The CTF fit has three parts. The power spectrum of the micrograph, a simulation of the power spectrum based on the fit, and an "astigmatism compensated" radial average of your data.

- 1. Compare the radial average to the power spectrum. If the bright and dark circles match, that means that the radial average represents our data well.
- 2. Compare the power spectrum and radial average to the simulated power spectrum. If they match, that means that the CTF parameters describe our data well and particles from the micrograph are likely to be useful for further analysis.

3. On manual-picking GUI:

Click CTF to examine the CTF fit Click pick to manually pick particles



Leicester Institute of Structural and Chemical Biology



#### Manual picking

This is the first place in the processing that we can bias our results. If we think we 'know' what our particle looks like, we will only pick particles (or noise) that looks like that. Instead, our job is to just pick things we think might be particles. Be as vague as possible!



- 1. Pick particles with the left mouse button
- Notice that the circles we get are more like dots. That's why we set 20 Å particle diameter
- 3. Middle-click in the center of the marker to remove picked particle

- 4. Right-click and Save STAR with coordinates on every micrograph before you close it, or you'll lose your picks!
- Now choose 500-1000 particles from micrographs with relatively high levels of defocus. Target 50-100 particles per micrograph.





#### **Particle sizing**

We need an initial guess for the size of our particles. If you have an idea from other sources (MAALS, PDB, etc,...) you can use that as an initial guess. But remember - your data is real, your initial guesses are not. You MUST choose your size based on the data!

In order to guess the size, we adjust the particle diameter parameter and *Continue!* The manual picking job, checking as we go.

000	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Display Colors Running
Import	Particle diameter (A): 100
CTF estimation	Scale for micrographs: 0.2
Manual picking	Sigma contrast: 3
Particle extraction	White value: 0
Subset selection	Black value: 0
3D initial model	Lowpass filter (A) 20
3D classification 3D auto-refine	Highpass filter (A) -1
3D multi-body	Pixel size (A) -1 7
Bayesian polishing Mask creation	Pick start-end coordinates helices? No
Join star files Particle subtraction	Scale for CTF image: 1

- 1. Click on your ManualPick/job005
- 2. Change the Particle diameter (A): 100
- 3. Click *Continue*!
- 4. Examine a few micrographs
- 5. Try changing the scaling to make the micrographs bigger





#### **Particle extraction**

Relion operates on particle images extracted from micrographs. We need to decide how big an image to extract, and by how much we want to reduce the resolution (Å/pixel). If your picking is not exactly correct (and it won't be,) choosing a bigger box size will let Relion optimize the position.

As a guideline (and assuming our protein is approximately spherical,) we want to select a box size twice as large as the largest dimension of the particle. We will then rescale it as small as we can (to increase the calculation speed,) whilst maintaining enough resolution for this stage of processing (3~4 Å/pix.)

	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O extract Helix Running
Import	micrograph STAR file: 1/ob003/micrographs_ctf.star P Browse
CTE estimation	
Manual picking	Input coordinates: pords_suffix_manualpick.star ? Browse
Auto-picking	
Particle extraction	OR re-extract refined particles? No 2
Subset selection	Refined particles STAR file: 7 Browse
2D classification	Reset the refined offsets to zero? No 2
3D classification	OB: re-center refined coordinates? No.
3D auto-refine	
3D multi-body	Recenter on - X, Y, Z (pix): 0 0 0 7
CTF refinement	Manually cot pixel cize? No
Bayesian polishing	
loin star files	Pixel size (A) 1
Particle subtraction	
Post-processing	
Local resolution	
External	
	Schedule Check command Kun!
I/O view Job actions	Current: Give_alias_here Display:

 On I/O tab: Micrograph STAR file: Select/all\_good\_micrographs /micrographs\_ctf.star

Input coordinates? ManualPick/jobXXX/ cords\_suffix\_manualpick.star

Here is a sample calculations for a 200 Å maximum dimension particle. We use a box size twice as big as our particle, and use at most 4 Å /pixel resolution. Note that Relion requires an even number of pixels:

$$\left(\frac{Particle Size Guess(\AA)}{Pixel Size(\AA/pix)}\right) \times 2 = BoxSize$$
$$\left(\frac{200 A}{0.885 \ \AA/pix}\right) \times 2 = 451.9 \approx 452 \ pix$$

$$BoxSize \times \left(\frac{Acquired \ Pixel \ Size \ (\mathring{A}/pix)}{Target \ Pixel \ Size \ (\mathring{A}/pix)}\right) = ExtractSize$$

$$452 \times \left(\frac{0.885 \ \mathring{A}/pix}{4 \ \mathring{A}/pix}\right) \approx 100 \ pix$$

At this stage, don't exceed 80 pixel boxes - even if it means operating at 10 Å/pixel. Our goal right now is to eliminate junk, which you can still do very easily at very low resolution. If you want to make your boxes smaller, you can go as low as 60 pixels in the first instance.



#### **Particle extraction**

	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O extract Helix Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification 3D initial model 3D auto-refine 3D auto-refine 3D auto-refine 3D auto-refine Bayesian polishing Mask creation Join star files Particle subtraction Particle subtraction Particle subtraction College Statement Distance Statement College Statement Distance Statement College Statement External	Particle box size (pix): 452 P Invert contrast? Yes Normalize particles? Yes Diameter background circle (pix): 200 P Studev for white dust removal: 1 P Studev for white dust removal: 1 P Rescale particles? Yes Re-scaled size (pixels): 100 P Schedule Check commant Run!

 On Extract tab: Particle box size: 452

Invert contrast? <mark>Yes</mark> (No for negative stain)

Normalize particles? <mark>Yes</mark>

Rescale particles? <mark>Yes</mark>

Re-scaled size (pixels): 64

0 0	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O extract Helix Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction 2D classification 3D initial model 3D classification 3D classification 4D classificati	Number of MPI procs: 3 P Submit to queue? Yes  P Queue submit commant: sput Auximum run time (d-h:m:s) 3-00:00:00 GPU count (0-4) 0 Standard submission script://net/prog/slurm/relion.sh 2 Browse Minimum dedicated cores per node: 1 Additional arguments: 2 P
EALEITIAI	Schedule Check command Run!
I/O view Job actions	Current: Give_alias_here Display:

On Running tab:
 Number of MPI processes:
 3 (16-32 for real data)

Submit to queue? <mark>Yes</mark>

Click Run! ONCE

Note: Relion will produce a warning for every micrograph missing particle location files. This is normal. On our cluster, this job takes ~10 sec.



#### **2D classification**

For 2D classification, more MPI processes are preferred over more threads because 2D classification requires less memory than 3D classification or refinement. Unlike Motion correction, CTF estimation or Auto-picking, classification and refinement processes need to be coordinated by a master process. Here, we use 1 GPU, 2 MPI processes per GPU and 2 threads per MPI process.

File lobs Schedules	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial	1. On I/O tab:
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection	Input images STAR file: Extract/job006/particles.star ) Browse Continue from here:	Input images STAR file: Extract/job006/particles.star
File Jobs Schedules File Jobs Schedules Impart Motion correction CTF estimation Auto-picking Particle extraction Subset selection D classification D ab classification D ab classification D ab classification File Jobs Schedules File Jobs Schedules Import Motion correction CTF estimation Manual picking Auto-picking D subset selection Subset selection D subset selection D auto-refine	RELION-3.1-beta-commit-7f8d1d: /home/lj/22/relion31-tutorial         IVO       CTF       Optimisation       Sampling       Helix       Compute       Running         Number of classes:       50       7       7       7       7       7         Number of leastes:       2       7       7       7       7       7         Number of treations:       25       7       7       7       7       7       7         Use fast subsets (for large data sets)?       No       2       7	<pre>2. On Optimization tab: Number of classes: 50 Mask diameter: 250 Limit resolution E-step to (A): 30 3. On Compute tab: Use parallel disc I/O: Yes</pre>
CTF refinement Bayesian polishing		Number of pooled particles: 30 Use GPU acceleration:
File Jobs Schedules	X: RELION-31-beta-commit-7/8d1d: /home/tjr/22//elion31-tutorial	Yes
Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection <b>2D classification</b> 3D initial model 3D classification 3D auto-refine	Number of MPI procs; 5 / 2 / 2 / 2 / 2 / 2 / 2 / 2 / 2 / 2 /	<ul> <li>4. On Running tab:</li> <li>Number of MPI procs:</li> <li>3 (see comment below)</li> </ul>
3D multi-body CTF refinement Bayesian polishing Mask creation Join star files Particle subtraction Post-processing	Standard submission script: [net/prog/slurm/relion.sh 2] Browse Minimum dedicated cores per node: 1 [] 2 Additional arguments: 2 2	Number of threads: 1 (see comment below) Submit to queue?
External	Schedule Check command Run!	Yes
		GPU count (0-4): <mark>1</mark>
		Click <mark>Run!</mark> ONCE

On our cluster, this job should take 2~3 mins.

If you encounter an out of memory error, you'll need to reduce the number of MPI per GPU. Increase the number of threads to keep (MPI × Threads)  $\leq$  (6 × GPUs). Always keep MPI × Threads  $\leq$  32.



### **2D classification analysis**

File Jobs Schedules	I/O CTF Optimisation Sampling Helix Compute Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification 3D initial model 3D classification 3D classificatio	Number of MPI procs: 5 1 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7
External I/O view Job actions	Schedule         Check command         Continue!           Current:         007: Class2D/job007/         Display:
	in: particles.star out: run_it025_data.star

un_1023_1110	Scale: 2	Min:	0	Max:	0
Sigma co	ontrast: 0	Color:	greysc	ale	\$
[	Display: rlnRe	eferencelr	nage		
Sort ima	ges on: rinCl	assDistrib	ution		
Reverse	sort? OAp	ply orienta	ations?	Read	whole stac
Nr. columns	: <mark>5</mark> Ori	scale: 1	Ma	k. nr. ima	ges: <mark>1000</mark>
				ſ	Display

- 1. Click on Display: drop-down
- 2. Choose out: run it025\_model.star

- 1. Change Scale to 2
- 2. Sort on rlnClassDistribution
- 3. Reverse sort
- Set the number of columns to 4
   Click Display!

		A clussedjoboooliun_tozo	jillouci.stui	
$\langle \phi \rangle$	٩	۲	$\langle \hat{\mathbf{S}} \rangle$	۲
0			0	$\bigcirc$
Ģ		۲	۲	۲

- Observations
  - · Central particles visible
  - Other, nearby particles seen on the edge
  - Some, but not many internal features visible in particles
- Interpretation ٠
  - Mask it too big alignment being driven by location of adjacent particles
- Next step ٠
  - Reduce mask size



#### **2D classification**

Rather than re-doing the 2D classification, we can change some parameters and continue the job. Because Mask diameter is one of the parameters we can change and continue the job, we'll do that by selecting the starting point, increasing the number of iterations by 5, and reducing the mask size.

	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O CTF Optimisation Sampling Helix Compute Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification 3D classification	Input images STAR file: [ <u>Extract/job006/particles.star</u> ] Browse Continue from here: <u>jo07/run_it025_optimiser.star</u> ] Browse

 Click on the previous 2D classification job in the completed jobs box to load the job

On I/O tab:
 Continue from here:
 (Select last iteration)



On Optimisation tab:
 Number of iterations:
 30

Mask Diameter: 160

Click Continue! ONCE

Once finished, Display: out: run ct 25 it030 model.star

	X C	lass2D/job008/run_ct25_it030_mc	odel.star	
۲	۲	۲	8	۲
۲	۲	۲	۲	8
۲	۲	۲	۲	

- Observations
  - Good classes
    - Single particle within mask
      - Particles entirely within mask
    - Internal features visible
    - Very dark background
  - Bad classes
    - 'Fuzzy'
    - · Density extends to mask
    - etc,...
  - Empty classes
    - Řelion tends to collapse 2D classes together. If this is a known issue. You can try turning on 'Ignore CTFs until first peak' on CTF tab.
- Interpretation
  - Good classes sufficient for reference-based particle picking
- Next step
  - · Reference-based picking



#### LoG picking

As an alternative to manual picking, we can use a generalized picker. These are less accurate than reference-based picking, but faster than manual picking. However, they do require accurate estimates of particle size to be effective, and the best way to do that is via manual picking and 2D classification.

Relion implements Laplacian-of-Gaussian (LoG) picking, which basically means it picks oval blobs within a set size range. Artificial-intelligence based pickers are more powerful, (and are available as External jobs now and likely to be included in Relion in the future,) but will use the same GUI approach.

The first step is to select 3-5 micrographs across the defocus range (using the OR select from micrographs.star option in the Subset selection job. Select from your CTF estimation job, display rlnCtfImage, and sort by rlnDefocusU. Give the job a sensible alias, like 3\_micrographs.)

	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Laplacian References autopicking Helix Running
Import Motion correction CTF estimation	Input micrographs for autopick: elect/job016/micrographs.star ? Browse Pixel size in micrographs (A) -1 7
Auto-picking Particle extraction Subset selection	2D references: 7 Browse OR: provide a 3D reference? No + 7
2D classification 3D initial model 3D classification	3D reference: 7 Browse Symmetry: C1
3D auto-refine 3D multi-body	3D angular sampling: 30 degrees + 7 OR: use Laplacian-of-Gaussian? Yes + 7
Bayesian polishing	

 File
 Jobs
 Schedules

 Import
 WO
 Laplacian
 References
 autopicking
 Helix
 Running

 Micin correction
 Min. diameter for LoG filter (A) 140
 0
 0
 0
 0

 Manual picking
 Max. diameter for LoG filter (A) 180
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0
 0

X RELION-3.1-beta-commit-7f8d1d: /home/tir22/relion31-tutorial File Jobs Schedules I/O Laplacian References autopicking Helix Running Import Motion correction CTF estimation Manual picking Particle extraction Subset selection 2D classification 3D initial model 3D classification 3D auto-refine Picking threshold: 0.05 Minimum inter-particle distance (A): 80 Maximum stddev noise: -1 Minimum avg noise: -999 ¢? Write FOM maps? Yes Read FOM maps? No 3D auto-refine 3D multi-body CTF refinement Shrink factor: 0 Use GPU acceleration? Yes - ? • ? D Bayesian polishing Mask creation Join star files Particle subtraction Which GPUs to use:

0 0	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Laplacian References autopicking Helix Running
Import Motion correction	Number of MPI procs: 3
CTF estimation	Submit to queue? Yes 🔷 ?
Auto-picking	Queue name: gpu ?
Particle extraction	Oueue submit command: sbatch
Subset selection	Maximum run time (d-h:m:s) 3-00:00:00
3D initial model	GPU count (0-4) 0
3D classification	Standard submission script. In strans follow failing ab.
3D auto-refine	Standard submission script. //net/prog/siurm/relion.sn / Browse
3D multi-body CTF refinement	Minimum dedicated cores per node: 1
Bayesian polishing	Additional arguments: ?
Mask creation	

1. On I/O tab:

Input micrographs for autopick: Select/3\_micrographs/micrographs.star

OR use Laplacian-of Gaussians? <mark>Yes</mark>

2. On Laplacian tab:

Min. diameter for LoG filter (A): 140 (a bit smaller that the minimum size of your particle)

Max. diameter for LoG filter (A): 180 (a bit bigger that the maximum size of your particle)

3. On autopicking tab: Write FOM maps? Yes Use GPU acceleration? No

4. On Running tab:
Number of MPI:
5 (or number of micrographs, whichever is fewer)
Submit to queue?

Yes

Click Run! ONCE



#### LoG picking

If you turn on Write FOM maps, Relion will write-out maps containing the Figure of Merit for each pixel (as generation of these maps is the most computationally expensive step, this is helpful for optimization of picking parameters.) If you then *Continue!* the job by reading the FOM maps, Relion with **replace** the particles in the previous round of picking with those chosen with the new parameters. Otherwise *Continue!* will ignore already picked micrographs and only pick new micrographs with the updated parameters.

We can now optimize our picking parameters.

	RELION ma	anual-pickin	g GUI		
Fil	2				
~	aphs/20170629_00028_frameImage.mrc	pick	349	CTF	10126.8
•	aphs/20170629_00029_frameImage.mrc	pick	478	CTF	8276.56
~	aphs/20170629_00049_frameImage.mrc	pick	362	CTF	13108.3
					1

Click on the finished job

 Display the results in the same was as for the 2D classification, choosing out: cords suffix autopick.star

2. Display all micrographs by clicking on pick Note: Relion takes the display parameters from the last Manual picking job. If you need to change something, select Manual picking, change the parameters you want, then choose the Jobs menu at the top and click Save job Settings



- Observations
  - Most (but not all,) of the particles have been picked
  - Some contamination picked
- Interpretation
  - Good enough for initial LoG pick. If not, see the section on reference-based picking to adjust the Default threshold.
- Next step
  - Extract particles and make 2D classes for references
     Please do this now, using a 320 Å box size (you need to convert to pixels) and 64 pixel re-scaled size.



#### Subset selection (Classes)

We now need to select classes to be our references for auto-picking our full dataset. Select the fewest classes  $(3\sim5)$  that have the <u>most varied viewpoints</u>. In the next step, they will be down-sampled to 20 Å, so resolution isn't important, shape is.

	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Class options Subsets Duplicates Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction <b>SUbset selection</b> 2D classification 3D classification 3D classification	Select classes from model.star: /job019/run_it025_model.star ? Browse OR select from micrographs.star: ? ? Browse OR select from particles.star: ? ? Browse OR select from picked coords: ? Browse

1.	On	I/O	tab:

Select classes from model.star: Class2D/job###/run\_it025\_model.star (Make sure all the other fields are empty.)

File Jobs Schedules	I/O Class options Subsets Duplicate	s Running	21
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection	Re-cent the class averages? Regroup the particles? Approximate nr of groups:	Yes No	5 ¢ 5 ¢ 5 7
2D classification 3D initial model 3D classification			

2. On Class options tab: Re-center the class averages? Yes



Display as for 2D classification



- 1. Select 3-5 classes:
  - Must look like protein projections
  - Try to get the most varied views
  - Don't just use the highest ranked
- 2. Right-click and Save selected classes
- 3. Close all display windows



#### **Reference-based picking**

The *de facto* standard way to pick your full dataset is via reference-based template matching. Now that we have good references (which are not always so easy to get,) we can use Relion's reference-based auto-picker and our 3-5 micrographs to tune the picking parameters.

00	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial	1. On I/O tab:
File Jobs Schedules	1/0 Laplacian References autopicking Helix Running	Input micrographs for autopick:
Import	Input micrographs for autopick: etc/iob016/micrographs star 7 Browse	input mitologiaphib for autopiton.
Motion correction CTF estimation	Pixel size in micrographs (A) -1	select/job###/micrographs.star
Manual picking		J J I
Particle extraction Subset selection	OR: provide a 3D reference? No	2D References:
2D classification 3D initial model 3D classification	3D reference: 2 Browse Symmetry: C1 2	Select/job###/class averages.star
3D auto-refine	3D angular sampling: 30 degrees 🔷 ?	
3D multi-body CTF refinement Pavesian polishing	OR: use Laplacian-of-Gaussian? No 🗘	OR use Laplacian-of Gaussians?
Mask creation		No

0 0	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Laplacian References autopicking Helix Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification	Lowpass filter references (A) 20 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2
3D classification 3D classification 3D auto-refine 3D multi-body CTF refinement Bayesian polishing Mask creation Join star files Particle subtraction	References have inverted contrast? Yes

2. On References tab: Pixel size in references (A):

(From extract job run.out file)

Mask Diameter (A): <mark>160</mark>

File Jobs Schedules	I/O Laplacian References autopicking Helix Running
Import	Picking threshold: 0.1
CTF estimation	Minimum inter-particle distance (A): 75
Manual picking	Maximum stddey poise: 1.2
Auto-picking Particle extraction	Minimum avri poise: -000
Subset selection	Mininan avg holse.
2D classification	Write FOM maps? No 🔷 ?
3D initial model	Read FOM maps? No 🗢 ?
3D auto-refine 3D multi-body	Shrink factor: 0 7
CTF refinement	Use GPU acceleration? Yes 🔷 ?
Bayesian polishing Mask creation	Which GPUs to use: 7
Join star files Particle subtraction	

3. On autopicking tab: Picking threshold 0.1

Minimum inter-particle distance (A): (half the estimated particle size)

Maximum stddev noise: -1 Write FOM maps? Yes

Shrink factor: 0 Use GPU acceleration?

X RELION-3.1-beta-commit-7f8d1d: /home/tir22/relion31-tutoria File Jobs Schedules I/O Laplacian References autopicking Helix Running Import Motion correction CTF estimation Manual picking Number of MPI procs: 1 - ? 10-Submit to queue? Yes \$ ? Queue name: gpu Queue submit command: sbatch Maximum run time (d-h:m:s) 3-00:00:00 GPU count (0-4) 1 Standard submission script: //net/prog/slurm/relion.sh Minimum dedicated cores per node: ? Browse Additional arguments: ?

4. On Running tab: Number of MPI: 1 (1 MPI/GPU) Submit to queue? Yes GPU count (0-4): 1

Click Run! ONCE

ta-commit-7f8d1d: /home/tjr22/relion31-tutorial References autopicking Helix Running Number of MPI procs: Submit to queue? Yes

Yes



## **Reference-based picking**

Click on the finished job, and display the results in the same was as for the LoG pick.





- Observations
  - All particles (except those on the edge) picked
  - Ice contamination picked
  - Some noise picked
- Interpretation
  - Good starting point for optimization
- Next steps
  - 1. Turn off Write FOM maps and turn on Read FOM Maps
  - 2. Turn **off** GPU acceleration and queue submission
  - 3. Increase Picking threshold
  - 4. Click Continue!
  - 5. Right-click on micrographs and reload coordinates
  - 6. Repeat steps 3-5 until all particles picked and minimal noise picked
- Observations
  - All particles (except those on the edge) picked
  - Ice contamination picked
  - Minimal noise picked
- Interpretation
  - Need to optimize background filters
- Next steps
  - 1. Repeat 3-6 above, varying Maximum stddev noise (start at 1.0)
  - Repeat 3-6 above varying Manimum average noise (start at -0.5)

Once you have optimized picking parameters, click on the Auto-picking job in the available jobs box to create a new job with the same parameters. Change the input micrographs to your all\_good\_micrographs job, turn off Read FOM Maps, turn on GPU acceleration, and pick your full dataset (if you get fewer than 10k particles, adjust your settings to increase the number of picks.) For real data, increase the number of MPI procs and GPUs to 2-4.



#### Particle set filtration via 2D classification

Now that we have finally done a full pick of our data, we can use 2D classification to filter out junk particles. Note that at this stage, our philosophy should be <u>to remove junk particles</u>. <u>NOT select good classes</u>. (note the numbers below are for example – yours may differ.)



6064 particles

- Extract your particles with the same box size as before, and 64 pixel re-scalled size (~30s for 3 MPI processes).
- Run 2D classification with 50 classes, 20 iterations and Use fast subsets (~10 mins).
- 3. Do Subset selection, turning on regrouping particles in the classes tab (50 groups.) Select anything that doesn't look awful.
- Run 2D classification on the particles from step 3 with 200 classes, 25 iterations. (~20 mins). With a "real" dataset, you only use fast-subsets when you have more than ~50k particles, so turn it off at this stage.
- 5. Re-select as above.
- Repeat 2D & Subset selection loop until no-more junk is evident. (5~20 mins per job). Here, one more iteration should be enough. In a full dataset you may need several more cycles.

Note that you can select a class, right-click and show metadata to find out how many particles are in each class. Classes with fewer than 100-200 particles will almost always look bad (fuzzy, blurred, lowresolution, etc,...) but may still be worth keeping.

If, at this stage, your classes are collapsing from hundreds to a dozen or so, try turning on Ignore CTFs until first peak? on the CTF tab. The output classes will look odd, but should separate junk better (this should be considered a "last resort," though.)



#### Initial model generation

In order to do 3D classification, we require an initial 3D model. Although there are many ways to do this, we'll only see one here: *ab-initio* model generation using *stochastic gradient descent*.

Note: you may be tempted to generate your initial model from (possibly hypothetical) PDB/atomic coordinate data. This is, by far, the easiest way to impose model bias on your data – and is strongly discouraged. We do not provide instructions on how to do this.

File lobs Schedules	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial	1. On I/O tab:
Import Motion correction CTF estimation Manual picking	Input images STAR file: Select/job031/particles.star ) Browse Continue from here:	<pre>Input images STAR file: Select/job###/particles.star</pre>
Manual picking		
File Jobs Schedules File Jobs Schedules Motion correction CTF estimation Manual picking Autopicking Particle extraction Subset selection Distributio	N RELION-3.1-beta-commit-7/8d1d: /home/tjr/22/relion31-tutorial         VO       CTF         Optimisation       SGD         Number of classes:       2         Mask diameter (A):       180         Flatten and enforce non-negative solvent?       2         Symmetry:       CL         Initial angular sampling:       15 degrees         Offset search range (pix):       6         Offset search step (pix):       2	2. On Optimization tab: Number of classes: 4 Mask diameter: 180 (A bit bigger than the mask for the 2D)
File Jobs Schedules Import Motion correction CTF estimation Mutopiking Particle extraction Subset selection 2D classification 3D auto-refine 3D multi-body GTF refinement Mask creation	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial         VO       CTF       Optimisation       SGD       Compute       Running         Use parallel disc VO?       Yes       \$ \$ \$       \$       \$ \$         Number of pooled particles:       30       • \$ \$       \$ \$       \$ \$         Skip padding?       Yes       \$ \$ \$ \$       \$ \$ \$ \$       \$ \$ \$ \$         Pre-read all particles into RAM?       No       \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$ \$	<pre>3. On Compute tab: Use parallel disc I/O?: Yes Number of pooled particles: 30 Skip padding? Yes Skip gridding? Yes Use GPU acceleration? Yes</pre>
File Jobs Schedules File Jobs Schedules Import Motion correction CTF estimation Manual picking Particle extraction 20 classification 30 autorefine 30 minut model 30 classification 30 autorefine 30 minut model 30 minut model	NELION-3.1-beta-commit-7/Bd1d: /home/tjr22/relion31-tutorial         VO       CTF       Optimisation       SGD       Compute       Running         Number of MPI procs:       3       0       0       0         Submit to queue?       Yes       0       0       0         Queue name:       gpu       0       0       0       0         Queue name:       istach       0<	<pre>4. On Running tab: Number of MPI procs: 3 (see comment below) Number of threads: 2 (see comment below) Submit to queue? Yes GPU count (0-4): 1</pre>
		Click <b>Run!</b> ONCE. This should take 20~60 mins.

3D processing (classification and refinement) requires more GPU memory than 2D processing, so we can use fewer MPI processes per GPU. For classification jobs, we can generally use 2 MPI/GPU (plus one master process.) Just like in 2D, though, if you encounter an out of memory error, you'll need to reduce the number of MPI per GPU. Increase the number of threads to keep (MPI × Threads)  $\leq$  (10 × GPUs). Always keep MPI × Threads  $\leq$  32.



#### **Initial model generation**

We will now assess the 4 initial models we generated using UCSF Chimera. Choosing the correct starting model is both important, and something of an art. In general, look for models that look like they may actually contain your 2D classes. They should also 'look' like protein, and have most of their density in the middle of the box, rather than as a shell on the outside of the box.



- 1. Download all four initial models from your processing directory to your local machine, and open UCSF Chimera
- Once Chimera loads, open the 4 models by File > Open, then browse to wherever you saved them, select all 4 run\_it300\_class00#.mrc files, then click open.
- 3. On the volume viewer window, do Features > Data Display Options, then change the Maximum number of histograms shown to 4
- 4. On the volume viewer window toolbar, click on data, then click on the dataset missing from the histograms
- 5. On the UCSF Chimera window toolbar, click Tools > Structure Comparison > Tile structures. Then click OK on the Tile Structures window
- On the UCSF Chimera window toolbar, click Favorites > Command line. In the Command line that appears at the bottom of the window, type in set independent and press enter.
- 7. Of the four Initial maps shown here, all four are reasonable, but the grey and light blue both look much better. are suitable and can be used. The grey map will be used going forward.



#### **3D classification**

Now that we have an initial model, we can finally start choosing 'good' particles, rather than eliminating junk. We do this via 3D classification.

	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
ile Jobs Schedules	I/O Reference CTF Optimisation Sampling Helix Compute Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification	Input images STAR file: Select/pb031/particles.star (?) Browse Continue from here: Reference map: b032/run_it300_class001.mrc (?) Browse Reference mask (optional): ?) Browse
2D initial model	
3D initial model	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
<ul> <li>D initial model</li> <li>O</li> <li>File Jobs Schedules</li> </ul>	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial

	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Reference CTF Optimisation Sampling Helix Compute Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification	Do CTF-correction? [Yes ♦]? Has reference been CTF-corrected? [Yes ♦]? Ignore CTFs until first peak? No ♦]?

	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Reference CTF Optimisation Sampling Helix Compute Running
Import Motion correction CTF estimation Manual picking	Number of classes: 4 1 Regularisation parameter T: 4 7
Auto-picking	Number of iterations: 25 25
Particle extraction Subset selection	Use fast subsets (for large data sets)? No
3D initial model	Mask diameter (A): 180 7
3D classification 3D auto-refine	Mask individual particles with zeros? Yes
3D multi-body CTF refinement	Limit resolution E-step to (A): -1
Bavesian polishing	

•	0		XR	ELION-3.1-bet	a-com	nit-7f8d1d: /hom	e/tjr22/relion	31-tutor	ial	
File	Jobs	Schedules	I/O	Reference	CTF	Optimisation	Sampling	Helix	Compute	Running
Imp	port				Us	e parallel disc I,	O? Yes		:	• 7
CTE	estima	tion		Nur	nber o	f pooled particl	es: 30 -	-0		- 2
Mar	nual picl	king				Skin paddir	m? Yes			• 0
Aut	o-pickin	g				Skip priddir	res ra? Vec			
Sub	ucie exu oset sele	raction		0		skip griddi	Ig: Tes			
2D	classific	ation		Pre-rea	id all p	articles into RA				
3D	initial m	odel		Copy parti	cles to	scratch directo	ry:			2
30	auto-ref	auon		Combine	e iterat	ions through di	sc? No			<del>\$</del> ?
3D	multi-bo	ody			llco	GPU acceleratio	m2 Voc			• •
CTF	refinen	nent			030	thick CBUs to u	in Tes			
Bay	/esian pi sk creati	olishing			v	mich GFUS to u	se:			
Loin	sk creat									
			N D	CLION 2.1 her		mit 760 d1d. /hon	e /tir22/relien	21 tuto	riel	
				LEIOIN-3.1-Dei	a-com	inic-710010. /iton	ie/grzz/renon	131-1410		
File	Jobs	Schedules	1/0	Reference	CTF	Optimisation	Sampling	Helix	Compute	Running
Im	port		<u> </u>		Nu	mber of MPI pr	(s) 2	<u> </u>		
Mo	tion cor	rection			Nu .	lumber of three		u m		38
Ma	r estima nual nic	kina			r	umber of threa	us: 4			
Aut	to-pickir	ng				Submit to que	ue? Yes			<b>÷</b> ?
Par	ticle ext	traction				Queue nan	ne: gpu			2
20	classific	ation			Queue	submit comma	nd: sbatch			?
3D	initial m	nodel		Мах	imum	run time (d-h:n	n:s) 3-00:00:	00		12
3D	classific	ation				GPU count (0	-4) 1			10
3D	auto-re	tine odv		Ct-	ndard	submission ser	int: Inot/prov	a /elu rao	relien ch	Drowso
СТІ	F refiner	nent		Sta	inuaru	submission scr	ipc. (net/prog	g/siurm,	renon.sh	( Browse
Bay	yesian p	olishing		Minimum c	learcat	ea cores per no	de: T			
Ma	sk creat	ion			Add	litional argume	nts:			2
Par	ticla cub	es								

1. On I/O tab: Input images STAR file: (last select job)

Reference Map: (your favorite 3D initial model)

2. On Reference tab:

Ref. map is on absolute greyscale? Yes (only for maps <u>directly</u> from Relion refinements, No for resized maps, etc,... If you're not sure, choose No)

3. On CTF tab: Has reference been CTF corrected? 1. Yes

4. On Optimization tab: Number of Classes: 4 Mask Diameter: 180

5. Compute and Running tabs: Same as for 3D initial model

Click Run! ONCE This job should take 15~20 mins.



#### Interlude: 2D classes vs 3D model

If you're not sure if your model may match the 2D classes, you can change the view in chimera and try to find the right perspective. Here's how class001 matches up



2D Class Chimera















- 1. In Chimera, select File > open and select the initial model.
- 2. On the Chimera window select Tools > Viewing Controls > Camera and change projection to orthographic.
- 3. Use the scroll wheel to zoom-out so you can see the entire map.
- 4. On the Volume Viewer, change the Style from 'surface' to 'solid'
- 5. Move the middle control point on the histogram to match the figure (left).
- 6. On the Relion window, go to your last 2D classification job, and Display: out: run\_it025\_model.star
- Set the Scale: to 2, Sort images on: rnlClassDistribution, and Reverse sort.
- 8. Click Display!
- Now you can use the left mouse button in Chimera to rotate the model around, comparing it to your 2D Classes from Relion



#### **3D classification**

To monitor our 3D classification, we look in three different places

1. Relion's output window:

Double-clicking on the window brings up the contents of the run.out file (no live-refresh, though)

Here we can see how the run time is progressing, the highest resolution map, the best angular accuracy achieved, and best offset accuracy achieved. Here's the output for the end of the run:

Auto-refine: Estimated accuracy angles= 1.759 degrees; offsets= 0.385 pixels CurrentResolution= 11 3786 Angstroms, which requires orientationSampling of at least 6.42857 degrees for a particle of diameter 200 Angstroms Oversampling= 0 NrHiddenVariableSamplingPoints= 387072 OrientationalSampling= 15 NrOrientations= 4608 TranslationalSampling= 2 NrTranslations= 21
Oversampling= 1 NrHiddenVariableSamplingPoints= 12386304 OrientationalSampling= 7.5 NrOrientations= 36864 TranslationalSampling= 1 NrTranslations= 84 ===================================
Expectation iteration 25 of 25 17/ 17 sec

Note that we've reached an alignment resolution (not model resolution) of 11.4 A. If we want higher resolution, we either need to use a smaller mask diameter (blue box) or finer orientational sampling (green box) – and hope our data support it.

2. Per-iteration model.star file

On the command line, typing

less Class3D/jobXXX/run\_it025\_model.star will show the details of iteration 25. Here is the most relevant section of that file:

data_model_classes					
<pre>loop_ rlnReferenceImage #1</pre>					
_rlnClassDistribution #2					
_rlnAccuracyRotations #3					
_rlnAccuracyTranslations #4					
_rlnEstimatedResolution #5					
_rlnOverallFourierCompleteness #6					
Class3D/job029/run_it025_class001.mrc	0.925191	1.759000	0.385000	11.378571	0.999540
Class3D/job029/run_it025_class002.mrc	0.010207	5.000000	1.246000	35.400000	0.826843
Class3D/job029/run_it025_class003.mrc	0.017992	3.595000	0.853000	26.550000	0.938008
Class3D/job029/run_it025_class004.mrc	0.046610	3.640000	0.890000	24.507692	0.988701

What we can see is that a single class has 92.5% of the particles, with an accuracy of 1.759 degrees and 0.358 pixels. The other classes are mostly empty, and any particles in those classes are most likely junk.

 We can also view the classes in UCSF Chimera (you should do that now) Along with information from the model.star file, you should be able to choose which class(es) to progress to refinement.

Note: This dataset typically doesn't require much 3D classification, but several rounds of Class3D/Select jobs are typically required.



#### **Detecting preferential orientation**

Preferential orientation, where the particles orient themselves on the grid, leads to anisotropic (not the same in all orientations,) resolution, and can make analysis of your map very difficult. Here's how to detect it.

- 1. Load a single model in chimera (here, we load class001.)
- 2. Using the File > load on the main Chimera window, open up the corresponding xxx\_angdist.bild file





The .bild files get rendered as columns around our model. The size and color of the column have the same meaning: the number of particles in that view. As you can see from this, there isn't a completely uniform distribution of views, however there are views from a number of different perspectives. This gives us confidence that there are enough views present and our data do not show preferential orientation.





#### Subsets in 3D

Selecting particles is slightly different between 2D and 3D.

- 1. We don't want to sort the classes (we choose which one(s) we want by inspecting the classes in Chimera.)
- 2. We need to ensure that we don't have any duplicate particles, which will interfere with 'gold-standard' resolution estimates later on

#### 6064 particles



5504 particles

1. Perform subset selection as you would for 2D Ensure that you don't sort the classes. Your results may look different – include all good classes.

2. Remove duplicates via another Subset selection job. On the I/O tab, choose the previous Select job in OR select from particles.star On the Class options tab, turn off re-group particles On the duplicates, turn on OR: remove duplicates?

Click Run!

Note: depending on the picking, it's possible that no particles are removed.



#### **Re-extraction**

Thus-far, the particles in this tutorial were extracted at 4.98 Å /pix, which gives a Nyquist resolution of 9.96 Å. Qualitatively, the 3D classifications appear to be giving maps at ~10Å, so we now want to reduce our Å/pix to continue. We'll start with a modest increase in resolution to 3.54Å/pix.

Classification and refinement jobs also improve the center position of our particles substantially. As a result, we no-longer choose our box size based on imprecise picking positions. Instead, we now choose the box size to capture the delocalized information due to the CTF of our microscope. Rosenthal and Henderson (2003) show what maximum useful size box we should use:

BoxSize in  $\mathring{A} = \frac{Maxmimum \ Defocus \ in \ \mathring{A}}{Constant \times Target \ Resolution \ in \ \mathring{A}} + Particle \ Size \ in \ \mathring{A}$ 

Rado Danev's guidelines for the value of Constant (25 for 300 keV, 20 for 200 keV, and 12.5 for 100 keV) are convenient here. (Get the defocus range by displaying the log file from the CTF estimation job):

BoxSize in 
$$\mathring{A} = \frac{13\ 000\mathring{A}}{20 \times 7.08\mathring{A}} + 160\mathring{A}$$

*BoxSize in*  $Å \approx 252 Å$ 

It's worth noting that larger boxes introduce more noise and higher defocus images have lower signal-to-noise at higher resolution, so this number should be considered an absolute maximum, and in our experience box sizes lower than this maximum will perform better once you start getting to Nyquist values below ~5Å. The minimum recommended box size is about 1.5 × particle size; we'll use a box size of 240Å.

	A RELIGIN-5.1-beta-commit-718010.7nome/gr22/religins1-totomai
ile Jobs Schedules	I/O extract Helix Running
Import	micrograph STAR file: elect/job004/micrographs.star 7 Browse
CTF estimation	
Manual picking	Input coordinates: Prowse
Auto-picking Particle extraction	OR re-extract refined particles? Yes \$
Subset selection	Refined particles STAR file: efine3D/job043/run_data.star ? Browse
2D classification 3D initial model	Reset the refined offsets to zero? No 2
3D classification	OR: re-center refined coordinates? Yes 2
3D auto-refine 3D multi-body	Recenter on - X, Y, Z (pix): 0 0 0 7
CTF refinement	Manually set pixel size?
Bayesian polishing Mask creation	Pixel size (A)
Join star files	FIXEI SIZE (A)
0.0	∑ RELION-31-beta-commit-7/8d1d: /home/tjr22/relion31-tutorial
ile Jobs Schedules	X RELION-31-beta-commit-7/8d1d: /home/tjr22/relion31-tutorial       VO     extract       Helix     Running
lle Jobs Schedules	X RELION-31-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial         I/O       extract         Helix       Running         Particle box size (pix):       328
ile jobs Schedules import Motion correction CTF estimation	X RELION-31-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial         MO       extract         Helix       Running         Particle box size (pix):       328         Invert contrast?       Yes         + 10
lie jobs Schedules import Motion correction CFF estimation Manual picking	X RELION-3.1-beta-commit-7/8d1d: /home/tjr22/relion31-tutorial         100       extract         Helix       Running         Particle box size (pix):       328         Invert contrast?       Yes         Normality and relative 2 Var.       100
ile Jobs Schedules Import Motion correction CFF estimation Manual picking Autopicking Particle extraction	X: RELION-31-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial         I/O       extract         Helix       Running         Particle box size (pix): [328       [8]         Invert contrast?       Yes       ? [7]         Normalize particles?       Yes       ? [7]
ile jobs Schedules import Motion correction CTF estimation Manual picking Autopicking Autopicking Di classification	X RELION-31-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial         VO       extract         Helix       Running         Particle box size (pix):       328         Invert contrast?       Yes         Yes       + ?         Normalize particles?       Yes         Diameter background circle (pix):       200         Yes       + ?         Yes </td

Rescale particles? Yes Re-scaled size (pixels): 164

3D classif

1. On I/O tab:

Micrograph STAR file: (Same as before) Input coordinates: (Empty) OR re-extract refined particles? Yes Refined particles STAR file: (Select particles.star file) OR: re-center refined coordinates? Yes

2. On Extract tab: Particle box size (pix): 272 Particle box size: 68



#### New extract parameters require new maps

We need maps to match our particle sizes. Relion 3.1 does allow us to mix box sizes and resolutions across datasets and initial models. However, it will use the box size and resolution of the first image/model is encounters, which for 3D classification or refinement is the initial model.

There are two ways to get higher resolution maps. The first way is to convert our current map using using relion\_image\_handler.

1. In the terminal window, find the Class3D output file:

ls Class3D/job###/run class001.mrc

2. Now call relion\_image\_handler, giving it the input (--i) and output (--o) file names, as well as the current (--angpix) and new resolution (--rescale\_angpix) and new box size (--new\_box) from the Extract job (note that the slashes at the ends of the line are for convenience. You can leave them out if your screen is wide enough). Here, I use an old pixel size of 3.186Å, a new pixel size of 3.54Å, and a new box size of 68 pixels. Make sure to use your own values here!

```
relion_image_handler \
--i Refine3D/jobXXX/run_class001.mrc \
--o Refine3D/jobXXX/run_class001_68box_3.54angpix.mrc \
--angpix 3.186 --rescale_angpix 3.54 --new_box 68
```

If you examine the particles datablock in the particles.star file from the extract job, you'll notice it has fields for rnlAnglePsi, rnlAngleRot, rnlAngleTilt, rnlOriginXAngst and rnlOriginYAngst, as well as rlnDefocusU, rnlDefocusV and rnlDefocusAngle – the three alignment angles, two shifts, and three defocus values required to do a reconstruction. So the other way to generate a new map is to do a reconstruction directly from the particles. We do this using relion\_reconstruct (making sure to include the --ctf flag to apply ctf correction.)

```
relion_reconstruct --ctf \
--i Extract/job027/particles.star \
--o Extract/job027/map.mrc
```

As a matter of routine, you should re-do the last classification job you did before reextraction. Often – the higher resolution will allow Relion to separate more junk. Please repeat the Class3D job using the particles and map from this extract job (when choosing the 3D map, you may have to go up two levels by clicking on the .../ twice to get to the actual Extract job), then do a Subset selection job before moving on.



#### **3D auto-refinement**

This is the first time we get to see how much resolution we may be able to get out of our data. The auto-refinement algorithm in Relion will progressively increase the precision of the alignment while avoiding overfitting. Unlike classification, refinement doesn't have a fixed number of iterations – it proceeds until it converges.

This is also the first time we're likely to run out of memory. If you see a ERROR: CudaCustomAllocator out of memory error, this means that the GPU has run out of RAM, and you have to change back to one MPI per GPU instead of the 2 MPI per GPU we've been using so-far.

Note: Our particles were extracted at 3.54 Å /pix, so our maximum resolution is 7.08 Å.

ila labs Schadulas	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial	1. On I/O tab:
lie jobs schedules	VO Reference CTF Optimisation Auto-sampling Helix Compute Running	
mport Motion correction	Input images STAR file: Select/job041/particles.star 7 Browse	input imges STAR file:
CTF estimation	Continue from here: 7 Browse	(3D select ioh)
Auto-picking	Reference map: b036/run_it025_class001.mrc ? Browse	
article extraction	Reference mask (optional):	Reference Map:
D classification		(vour 3D class)
		(your 3D class)
	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial	0. On Defense take
ile Jobs Schedules	I/O Reference CTF Optimisation Auto-sampling Helix Compute Running	2. On Reference tab:
Import	Do CTF-correction? Yes \$ ?	Ref. map is on absolute grevscale?
CTF estimation	Has reference been CTF-corrected? Yes 🗘 ?	
Manual picking Auto-picking	Ignore CTFs until first peak? No	Yes (only for maps directly from Relion
Particle extraction		refinements No for resized mans etc.
Subset selection		rennements, no for resized maps, etc,)
	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial	2 On CTE tab
ile Jobs Schedules	I/O Reference CTF Optimisation Auto-sampling Helix Compute Running	
nport Action correction	Mask diameter (A): 200 7	Has reference been CTF corrected?
CTF estimation	Mask individual particles with zeros? Yes	
Manual picking Auto-picking	Lise solvent-flattened FSCs? No.	1. <mark>165</mark>
Particle extraction		
ile lobs Schedules	X RELION-3.1-beta-commit-7/8d1d: /home/tjr22/relion31-tutorial	4 On Optimization tab:
Import		4. On Optimization tab.
Motion correction	Use parallel disc I/O? Yes	Mask Diameter:
CTF estimation Manual picking	Number of pooled particles: 30 - 7	
Auto-picking	Skip padding? Yes 7	180
Subset selection	Pre-read all particles into RAM? No 2 7	
2D classification 3D initial model	Copy particles to scratch directory:	
3D classification	Combine iterations through disc? No	
3D multi-body	Use GPU acceleration? Yes	5 On Compute teh: Same on Class?D
CTF refinement Bayesian polishing	Which GPUs to use:	5. On Compute tab. Same as Class3D
Mask creation Join star files	0	
	RELION-3.1-beta-commit-7f8d1d: /home/tir/22/relion31-tutorial	
le Jobs Schedules	I/O Reference CTF Optimisation Auto-sampling Helix Compute Running	6. On Running tab: Same as Class3D
mport		
Motion correction	Number of MPI procs: 3	Note: The description of running parameters for
Motion correction CTF estimation	Number of MPI procs:         3         1         2 <th2< th=""> <th2< th="">         2         <th2< th=""></th2<></th2<></th2<>	Note: The description of running parameters for
Motion correction CTF estimation Manual picking Auto-picking	Number of MPI procs: 3 1 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2	Note: The description of running parameters for 3D classification applies here as well.
Motion correction CTF estimation Manual picking Auto-picking Particle extraction	Number of MPI procs: 3 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7	Note: The description of running parameters for 3D classification applies here as well.
Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification	Number of MPI procs: 3 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7	Note: The description of running parameters for 3D classification applies here as well.
Motion correction CTF estimation Manual picking Particle extraction Subset selection 2D classification 3D initial model 3D classification	Number of MPI procs: 3 7 Number of threads: 2 7 Submit to queue? Yes 7 Queue name: gpu 7 Queue submit command: sbatch 7 Maximum run time (d-h:m:s) 3-00:00:00 7	Note: The description of running parameters for 3D classification applies here as well.
Motion correction CTF estimation Manual picking Auto-picking article extraction Subset selection 20 classification 4D initial model 4D classification 4D autorefine	Number of MPI procs: 3 7 7 Number of threads: 2 7 7 Submit to queue? Yes \$ 7 Queue name: gpu 7 Queue submit command: sbatch 7 Maximum run time (d-hrms) 3-00:00:00 7 GPU count (0-4) 1 7	Note: The description of running parameters for 3D classification applies here as well.
Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification 3D initial model 3D classification 3D autorefine 3D multi-body CTF refinement	Number of MPI procs: 3 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7 7	Note: The description of running parameters for 3D classification applies here as well.
Motion correction CTF estimation Manual picking Auto-picking Auto-picking Particle extraction SUbset selection 2D classification 3D initial model 3D classification 3D autorefine 3D autorefine Bayesing polishing State and State State State State State State State State State CTF refinement Bayesing polishing	Number of MPI procs: 3 7 7 Number of MPI procs: 3 7 Submit to queue? Yes 7 Queue name: 9pu 7 Queue submit command: sbatch 7 Maximum run time (6-h:ms:) 3-00:00:00 7 GPU count (0-4) 1 7 Standard submission script: //net/prog/slum//relion.sh 7 Minimum dedicated cores per node: 1 7 7	Note: The description of running parameters for 3D classification applies here as well.
Motion correction CTF estimation Manual picking Auto-picking Auto-picking Particle extraction Subset selection 2D classification 3D initial model 3D classification 2D autorcline Bayesian polishing Mask creation Join star files	Number of MPI procs: 3 7 7 Number of MPI procs: 3 7 7 Submit to queue? Yes ? 7 Queue name: ppu ? 7 Queue submit command: sbatch ? 7 Maximum run time (d-hrm:s) 3-00:00:00 ? 7 GPU count (0-4) 1 ? 7 Standard submission script: //netprog/slurm/relion.sh ? Browse Minimum dedicated cores per node: 1 7 Additional arguments: 7	Note: The description of running parameters for 3D classification applies here as well. Click <b>Run!</b> ONCE
Motion correction CTF estimation Manual picking Autopicking Particle extraction 2D binkial model DD binkial model DD binkial model DD binkial model DD binkial model DD binkial DD binkial	Number of MPI procs: 3 7 7 Number of MPI procs: 3 7 Submit to quee? Yes 7 Queue name: gpu 7 Queue submit command: sbatch 7 Maximum run time (4/h.ms): 300:00:00 GPU count (0-4) 1 7 Standard submission script: //net/prog/slum/relion.sh 7 Minimum dedicated cores per node: 1 7 Additional arguments: 7	Note: The description of running parameters for 3D classification applies here as well. Click Run! ONCE

When the job finishes, it will report the unmasked resolution (i.e. the resolution of the entire volume. To find out how well we really did, we need to post-process our map (which means we need a mask.)



#### **Mask generation**

Generation of masks can be done many ways. The most common is via UCSF Chimera. Here, we generate a solvent mask from our auto-refined map to use in post-processing.

Mask generation proceeds in 4 steps.

- 1. The map is low-pass filtered to remove any high-resolution features
- 2. It is then binarised by setting all values below a threshold to zero, and all those above to one
- 3. The map is expanded by some number of pixels
- 4. A  $[0 \pi/2]$  cosine-shaped density is added to soften the edge of the map

In order to properly determine the threshold for binarization, it is helpful to first low-pass filter your map before downloading it and loading in into UCSF Chimera:

```
relion_image_handler \
--i Refine3D/job###/run_class001.mrc \
--o Refine3D/job###/run_class001_15Alowpass.mrc \
--lowpass 15
```

In order to properly determine the threshold for binarization, it is helpful to first low-pass filter your map before loading in into UCSF Chimera.

Now load run\_class001.mrc and run\_class001\_15Alowpass.mrc into Chimera. Set the step sizes for both maps to 1. Use the unfiltered map to see your density and set the unfiltered map level as high as possible before seeing 'dust' (noise) come up.

#### 0.05: Too high



File Jobs Schedules

File Jobs Schedules

Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification 3D initial model 3D classification

Import Motion correction CTF estimation Manual picking

#### 0.02: Just Right



0.015: Too low



1. On I/O tab: Input 3D map: Refine3D/job###/run\_class001.mrc

2. On Mask tab: Initial binarization threshold: (Just Right value)

Extend binary map this many pixels: 2

Add a soft edge of this many pixels: 5-10

3. On Running tab: Number of threads:

(Or you can submit to the cluster using up to 16 threads



X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial

X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial

Extend binary map this many pixels: 2 Add a soft-edge of this many pixels: 10

Lowpass filter map (A) 15

Pixel size (A) -1

Input 3D map: e3D/job043/run\_class001.mrc ? Browse

I/O Mask Helix Running

I/O Mask Helix Running



#### **Post-processing**

Relion can auto-sharpen our map, and report the resolution within the solvent mask we just generated, using the Post-process job. Note that this only generates reliable sharpening if your resolution is significantly better than 10 Å (which, at 6.8 Å, ours is already.)

	🔀 RELION-3.1-beta-commit-f2c3d8: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Sharpen Filter Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification 2D initial model	One of the 2 unfiltered half-maps: run_half1_class001_unfil.mrc (?) Browse Solvent mask: rtaskCreate/job044/mask.mrc (?) Browse Calibrated pixel size (A) (3.2037

1. On I/O tab:

One of the 2 unfiltered half-maps: Refine3D/job###/run\_half1\_class001\_unfil.mrc

Solvent mask: MaskCreate/job###/mask.mrc

	X RELION-3.1-beta-commit-f2c3d8: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O Sharpen Filter Running
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection 2D classification 3D clintial model 3D clintial model 3D auto-refine	Estimate B-factor automatically? Yes   Estimate B-factor automatically? Yes   Lowest resolution for auto-8 fit (A): 10   Use provided B-factor? No   MTF of the detector (STAR file) Imon/MTF/mtf k2_300kV.star ? Browse  Original detector pixel size: 0.885   ?

2. On Sharpen tab:

MTF of the detector (STAR file) /net/common/MTF/mtf k2 300kV.star

Original detector pixel size: 0.885

You can choose to submit to the queue or not. Click *Run!* ONCE

The output of this job has lots of information. Most notably, the b-factor Relion applies for sharpening is -229, and the calculated final resolution is 7.08 Å, which is exactly Nyquist.

As with many other Relion jobs, Post-processing generates a PDF output file. In this case, various FSC (resolution) and Guinier (sharpening) curves are plotted. *Please* Display *out: logfile.pdf* now.



The most important line in the FSC plot is the black line, show the FSC for the masked map after correction for masking effects. Here, we can see that the FSC doesn't reach zero, indicating that our data will probably support higher resolutions than the current Nyquist resolution.



The most important line in the Guinier plot is the red line. Ideally we want it to remain flat at higher resolutions, before falling off at the highest resolution. Here, you can see how there's a small rise at the higher resolutions, indicating very slightly over sharpening (this degree of over-sharpening is not significant.)



#### Next steps

From the output of the latest 3D auto-refine and Post-process jobs, we can see that the map is at (or very close to) the Nyquist resolution, and the FSC didn't fall to zero. We should therefore be able to extract more resolution from our data.

- 1. Re-extract the particles from the run\_data.star file from the Refine3D job. Remember that Relion requires even box sizes, so increase the Particle box size to 280 and the Re-scaled size (pixels) to 140 pixels.
- 2. Reconstruct a map from these particles.
- 3. Skip the Class3D job for this tutorial. Run a new Refine3D job with the newly extracted particles and map. (Note how much longer this job takes, even with only ~5k particles!)
- 4. Make a mask and Post-process the new map.

This time, the b-factor is -59 and the map went to 3.9Å. The log file shows that the FSC went to zero – so going to a finer Å/pix will not deliver more resolution (and may make it worse as more noise than signal will be introduced).



If you examine the penultimate iteration in the Refine3D job, you should notice that the estimated resolution before the two half-maps were combined was 8~9Å, which then dropped to 4~5Å once the maps were combined. This is strongly suggestive that signal-to-noise is currently limiting our resolution, rather than alignment quality.

So it seems clear that we need to either remove more bad particles (Using 2D or 3D classification, with or without <u>image alignment</u> (look in the <u>Sampling</u> tab); or we need to add more good particles.

In a real dataset, the next step would be to filter out more bad particles, but we'll skip that now. DO NOT skip this step on real data - filtering out bad particles is <u>much more powerful</u> than adding new good particles!

There are two ways of adding good particles: Collecting more data or use symmetry.



### **Using symmetry**

The main problem with symmetry is that is assumes your particle is symmetric. Any nonsymmetry in your particles will reduce the resolution of your map. Applying symmetry will also cause you to lose any biologically relevant asymmetry information. Although using symmetry often gets you a prettier map, it may also be the wrong thing to do. You have been warned!

There are two steps required to use symmetry (see note below before actually doing these steps):

- 1. Align your map to the appropriate symmetry axis (axes).
  - 1. In the terminal, call relion\_align\_summetry, giving it the input (--i) and output (--o) file names, as well as the symmetry operator (--sym), in this case, D2.

```
relion_align_symmetry \
```

```
--i Refine3D/job###/run_class001.mrc \
```

--o Refine3D/job###/run\_class001\_D2-aligned.mrc \

--sym D2

2. Although strictly unnecessary, you can then apply the symmetry to your map using relion\_image\_handler:

```
relion\_image\_handler \ \backslash
```

- --i Refine3D/job###/run\_class001\_D2-aligned.mrc \
- --o Refine3D/job###/run\_class001\_D2-aligned\_D2-applied.mrc \

--sym D2

It may be interesting to apply symmetry to you unaligned map as well. *Please do that now and examine the results in Chimera.* 

 Repeat your Refinements with the new symmetry-aligned starting map and forcing symmetry in 3D auto-refine (symmetry in the Reference tab of the 3D auto-refine job)

The tutorial data would now go to Nyquist resolution (3.54Å.) Normally, you would apply symmetry, check the resolution of the auto-refinement, and then re-extract the particles and re-scale the maps. For expediency, please do the re-extraction, rescaling, symmetrical refinement and post-processing now (use a 280 pixel box and down-scale particles to 200 pixels.)

You should now have a map at 3.3~3.4Å resolution (well above the new Nyquist resolution of ~2.5Å. In a real dataset, now would be a good time to remove more bad particles. For this tutorial, we will proceed directly to refinement of the CTF parameters and per-particle motion correction. Because it's generally quicker, we start with a basic CTF refinement.



#### **CTF refinement**

CTF refinement works by comparing your particles to your map – so the better the map, the better CTF parameter estimates you can generate – which should get you a better map. As a result, you may need to do 2 or more cycles of CTF refinement + 3D auto-refine.

At low resolutions (worse that ~4.5 A,) only the defocus and large values of astigmatism have a significant impact on your map. Once you start to go to higher resolutions, other optical aberrations of the microscope can become limiting. Relion has the ability to estimate these additional aberrations using the CTF refinement job. Subsequent rounds of refinement will then take these estimates into account.

The first step is to refine the defocus estimates for all of our particles. In addition, we can estimate the beamtilt and trefoil (3<sup>rd</sup> order) aberrations. Other aberration estimations are possible, but we don't recommend using them at this resolution.

File Jobs Schedules	X       RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial         VO       Fit       Running         Particles (from Refine3D): [Polish/job052/shiny.star	<ol> <li>On I/O tab: Particles (from Refine3D): (use the run_data.star file from your refine job) Postprocess STAR file: (matching post-process job from refine job)</li> </ol>
File Jobs Schedules File Jobs Schedules Import Motion correction CTF estimation Manual picking Particle extraction Subset selection 20 classification 3D auto-refine 3D classification 3D auto-refine 3D auto-refi	RELION-3.1-beta-commit-7f8d1d: /home/ljr22/relion31-tutorial         VO       Fit         Bunning         Estimate (anisotropic) magnification?       No         Perform CTF parameter fitting?       Yes         Fit astigmatism?       Per-micrograph         Fit astigmatism?       No         Fit bef-actor?       No         Fit bef-actor?       No         Fit bef-actor?       No         Estimate beamtit?       Yes         Also estimate trefoil?       Yes         Estimate 4th order aberrations?       Yes         Minimum resolution for fits (A):       30	<pre>2. On Fit tab: Estimate (anisotropic magnification? No Perform CTF parameter fitting? Yes Fit Defocus? Per-particle Fit astigmatism? Per-micrograph Estimate beamtilt? Yes Also estimate trefoil? Yes</pre>
File Jobs Schedules File Jobs Schedules Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction 2D classification 3D classification Baselian polishing Mask creation Join Starf files	N RELION-3.1-beta-commit-7f8d1d: /home/ljr22/relion31-tutorial         VO       Fit         Number of MPI procs:       1         Number of threads:       2         Submit to queue?       2         Queue name:       2         Queue submit command:       2         Maximum run time (d-hzms) 3:00:000       2         GPU count (0:4)       2         Standard submission script:       1         Minimum dedicated cores per node:       2         Minimum dedicated cores per node:       2         Maximum run time:       2	<pre>3. On Running tab: Number of MPI procs: 1 (4 for real data) Number of Threads: 5 Submit to queue? Yes Click Run ( ONCE This should take 1~2 mins)</pre>

If your map is low-resolution, especially noisy, or your particles are especially small, fitting the beamtilt and trefoil may actually decrease the resolution of your reconstruction. Because polishing and CTF refinement are iterative, it is often best to fit just the defocus and astigmatism on the first iteration and only fit the others once the map no-longer improves. As phase-shift only makes sense for phase-plate data, and even then only plays a major role at low resolutions, so it's common not to fit the phase-shift at all.

Now repeat the 3D auto-refine and post-process jobs to see if there's any improvement (make sure to use the CtfRefine job as the input particles!)



### **CTF refinement (output)**

The logfile from the CTF refinement job contains a wealth of information.

#### **Defocus refinement**

There is an even distribution from Blue to Yellow in most micrographs, indicating that the particles are distributed across the thickness of the ice. Micrographs with particles localized to the ice/water interface show strongly blue and strongly yellow particles with nothing between.



Several micrographs were tilted. Note that the particles on, for example micrograph 31, go from mostly blue to mostly gold

#### <u>Asymmetrical Aberrations</u> Beam tilt is evident from the lower scale asymmetric plot (top left plot) as a blue-red transition across the plot. No trefoil aberrations are evident (they would be triangular blue-red shapes.)



The fitted Zernike polynomials (below) used for correction seem to match the measured aberrations (above).



The fitted beam tilt is 0.03 x -1.17 mRad. (found in the particles\_ctf\_refine.star file.)

#### Symmetrical Aberrations

A small Cs error is evident from the lower scale symmetric plot (top left plot) as a blue ring at higher resolutions (inner to outer.) No tetrafoil aberrations are evident (they would be square blue-red shapes.)



The fitted Zernike polynomials (below) used for correction seem to match the measured aberrations (above).





## **Anisotropic Magnification (CTF refinement)**

We can also estimate the anisotropic magnification using the CTF refinement job. Combining the other estimation tasks with this one is not possible, so we simply take the output of the previous job and use it as input here. (We actually run this later in the tutorial.)

000	X RELION-3.1-beta-commit-7f8d1d: /hom	e/tjr22/relion31-tutorial		
File Jobs Schedules	I/O Fit Running			
Import Motion correction CTF estimation Manual picking	Particles (from Refine3D): <u>b053/particles_ctf_refine.star</u> ] Postprocess STAR file:[ <u>cess/job051/postprocess.star</u> ] Browse			
	K RELION-3.1-beta-commit-7f8d1d: /home/tjr	22/relion31-tutorial		
ile Jobs Schedules	I/O Fit Running			
Import Motion correction	Estimate (anisotropic) magnification?	Yes 🗘	)?	
CTF estimation Manual picking	Perform CTF parameter fitting?	No	2	
Auto-picking	Fit defocus?	Per-particle \$	2	
Particle extraction	Fit astigmatism?	Per-micrograph 🔶	2	
2D classification	Fit B-factor?	No	2	
3D initial model	Fit phase-shift?	No 🗘	2	
3D auto-refine				
3D multi-body	Estimate beamtilt?	No		
CTF refinement Bayesian polishing	Also estimate trefoil?	Yes 🌲	2	
Aask creation	Estimate 4th order aberrations?	No 🌢	2	
oin star files		¥.		
Particle subtraction Post-processing	Minimum resolution for fits (A):	30	2	

 On I/O tab: Particles (from Refine3D): (use the file from your previous ctf refinement job) Postprocess STAR file: (same as for previous polishing job)

 On Fit tab: Turn everything off, then:

Estimate (anisotropic magnification? Yes

Click Run! ONCE. This should take 1~2 mins.

Once the job finishes, examine the resulting logfile. Similar to the aberrations in the previous step, we have measured X and Y distortion in the top plots, and fitted corrections in the bottom plots.



Although it is clear from the logfile that the CTF refinement job can detect anisotropic magnification in the data, the output matrix looks (something) like this:

 $\begin{bmatrix} 1.001341 & 0.001665 \\ 0.001609 & 0.997878 \end{bmatrix}$ 

With anisotropic magnification on the top left and bottom right, and skew on the top right and bottom left. This indicates that there's only  $\sim 0.3\%$  anisotropic magnification or skew, which is insignificant at this particle size.



#### **Bayesian polishing**

It is convenient to think of Bayesian polishing as a re-extraction directly from the original micrograph movies. Unlike standard motion correction, polishing fits map projection to each frame to estimate the movement (as well as the per-frame resolution.) So, just like CTF refinement, the better the map, the better the polishing – which should give you a better map.

While it requires you to tune several parameters for your particular dataset for **best** performance, you can do at least one round of polishing using the default parameters before you have to do the fitting (which takes a long time!)

Determining precise extraction parameters is an art. For now, we'll aim for 3Å target resolution – suggesting we should use a 426 pixel box down-sampled to 252 pixels. But see the next slide on how this is done in real projects.

File Jobs Schedules File Jobs Schedules Import Motion correction CTF estimation Manual picking Autopicking Particle extraction Subset selection Distail model Distail model Dis	N         RELION-31-beta-commit-7f8dfd: /home/tjr22/relion31-tutorial           VO         Train         Polish         Running           Micrographs (from MotionCorr):         [2/corrected_micrographs.star]         [2] Browse           Particles (from Refine3D or CtfRefine):         efine3D/job049/run_data.star]         [2] Browse           Postprocess STAR file:         cess/job051/postprocess.star]         [2] Browse           First movie frame:         1         7           Last movie frame:         1         7           Re-scaled size (pixels):         1         7	<pre>1. On I/O tab: (Fill in appropriate values) Extraction size (pix in unbinned movie): 426 Re-scaled size (pixels): 252</pre>
		2. On Train tab: Train optimal parameters? <mark>No</mark>
File Jobs Schedules File Jobs Schedules Import Maton Sometion Manual picking Auto-picking Particle extraction Subset selection 20 classification 30 initial model 30 basic-ordine 30 basic-ord	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial         VO       Train         Polish       Running         Perform particle polishing?       ?         Optimised parameter file:       ?         OR use your own parameters?       ?         Sigma for velocity (A/dose):       ?         Sigma for acceleration (A/dose):       ?         Sigma for acceleration (A/dose):       ?         Minimum resolution for B-factor fit (A):       1         Maximum resolution for B-factor fit (A):       1	<pre>3. On Polish tab: Perform particle polishing? Yes OR use your own parameters? Yes Submit to queue? Yes</pre>
File Jobs Schedules File Jobs Schedules Import Motion correction CTF estimation Manual picking Particle extraction Subset selection 3D initial model 3D classification 3D auto-refine 3D auto-refine 3D auto-refine 3D auto-refine 3D auto-refine 3D auto-refine Auto-refine Mask creation Particle subtraction	Image: Second State Sta	<pre>4. On Running tab: Number of MPI procs: 1 (4 for real data, or fewer if you get a memory error) Number of Threads: 5 Submit to queue? Yes For real data: Queue submit command: sbatchexclusive Click Run! ONCE. This should take ~10 mins.</pre>

Now repeat the 3D auto-refine and post-process jobs to see if there's any improvement (make sure to reconstruct a new map and align it to D2, and use the Polishing job shiny.mrc file as the input particles!)

Now that we're trying to squeeze the best out of our data, we want to set Mask individual particles with zeros? to No in the 3D auto-refine Optimization tab and set the Skip Padding? Option to No in the Compute tab.



#### Choosing your box size

This tutorial has so-far presented a fairly complex, math-heavy way of choosing your box size. In reality, most people just choose a box size twice the size of their particle. For example, for have a 160 Å particle at 1.5 Å/pix, a typical final box size would be ~214 pixels, with most users choosing 200, 220, or sometimes 256 pixels.

Ultimately, the box size you choose for your final, highest-resolution, map(s) is a compromise between two factors:

- · Bigger boxes incorporate more delocalization information
- Bigger boxes incorporate more noise.

... so the only real option is to try several box sizes out and see which gives the best map (or, if you're looking for a quick proxy, the least-negative fitted sharpening b-factor.) Relion does, however, provide one hint as to how big our box sizes should be.

- 1. Plot the delocalization for some particular particle set, including the radius of your particle (in this case 160Å)
  - relion\_plot\_delocalisation \
    --i Polish/job###/shiny.star \

```
--o Polish/job###/deloc \
```

```
--rad 160
```

 Open the plot using evince (or download to view with your pdf viewer of choise) evince Polish/job###/deloc.eps



The black line shows how many pixels of information you exclude for the center of your particle at some particular box size. The blue line is a bit more complicated. It shows the 'worst case'. You can think of it like this: for a spherical particle, how many of the particle images have all the delocalized information captured before you get to the nearest edge.

So, for our box size of 252 pixels, we have all of the information from the center of our particles, but only ~15% of the images have all of the highest resolution information from the edge of the particle. This may sound bad, but recall that the delocalization is circular, so most of the highest-resolution information will still be captured in most of the particle images, along with virtually all of the lower-resolution information.

The take-home message from the plot is this: There is virtually no information for boxes over 300 pixels, and none at all for boxes greater than 350.



#### Next steps

After post-processing, the map is at 3.09Å with a fitted b-factor sharpening of -22.71.

Now we repeat the CtfRefine – Polish – 3D auto-refine cycle until we get no improvement.

- 1. CtfRefine
  - Turn on 4<sup>th</sup> order aberration estimation normally we wouldn't do this until we were below 3Å.
- 2. CtfRefine (anisotropic magnification)
  - Use particles from previous CtfRefine job
  - Turn on Estimate (anisotropic) magnification. This will turn off all the other estimations.
- 3. Run a 3D auto-refine job (optional)
- 4. Run a Bayesian Polishing job using default settings. For a 2.6Å target resolution, we use a box size of 464 pixels down-sampled to 316 pixels.
- 5. Run a 3D auto-refine job and post-process. You should now have a map at <3.0Å and fitted b-factor sharpening of > -20 (i.e. -19)
- Repeat the above cycle, using the new map.
  - Train the Polishing job (this takes a long time.)
    - Normally, you would use 5-15K particles for training. The more particles, the better the parameter estimates but you may run out of memory. If this happens, reduce the number of particles and try again. If this still doesn't work, set --sbs in the Additional arguments field in the Running tab (but this will make training even slower!)
    - In this case, because we only have 4-5k particles, use 3000 particles.
    - On the Running tab, set:
      - Queue submit command to sbatch --exclusive -C skylake to ensure you use the new nodes with more memory and you reserve all the memory on the node.
      - MPI to 1 and Threads to 16.
  - Use the trained parameters by pointing the Optimised parameter file: to your training job output.
- Once you've done the post-processing, repeat the Refine3D and post-processing jobs:
  - Reference map: run\_class001.mrc from previous Refine3D job
  - · Reference mask: soft mask derived from the reference map
- Repeat the above, turning on solvent-flattened FSCs on the Optimization tab of the Refine3D job.
- The resolution (sharpening b-factor) obtained in writing this tutorial, after post-processing, were:
  - Refine3D: 2.85Å (-19.84)
  - Refine3D with masking: 2.85Å (-18.25)
  - Refine3D with masking and solvent-flattened FSCs: 2.85Å (-18.37)
- People often want their final map to be a 'full' resolution map. If you try this, you'll find a box size of 464 pixels is too big for the GPU memory. Instead, repeat the final Polish, Refine and Post-processing jobs with a box size of pix (and no down-sampling).



#### **Local Resolution**

Once you have obtained a density map, it can be useful to display the local resolution estimate to get a better feel for your results. Here, we use Relion's built-in local resolution estimator for the calculation, then UCSF Chimera for display.

	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial		
File Jobs Schedules	I/O ResMap Relion Running		
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction Subset selection	One of the 2 unfiltered half-maps: [run_half1_class001_unfil.mrc] [2] Browse User-provided solvent mask: [4askCreaterjob071/mask.mrc] [7] Browse Calibrated pixel size (A) [1_2744 []]		

00	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O ResMap Relion Running
Import Motion correction CTF estimation Manual nicking	Use ResMap? No + ? ResMap executable: esMap/ResMap-11.4-ilnux64 ? Browse

000	RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial	
File Jobs Schedules	I/O ResMap Relion Running	
Import Motion correction CTF estimation Manual picking Auto-picking Particle extraction	Use Relion? Yes + ? User-provided B-factor: -26.7 - 2 MTF of the detector (STAR file) ? Browse	

0 0	X RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial
File Jobs Schedules	I/O ResMap Relion Running
Import	Number of MPI procs: 5
CTF estimation	Number of threads: 2 - 7
Manual picking	Submit to queue? Yes \$
Particle extraction	Queue name: gpu ?
2D classification	Queue submit command: sbatch
3D initial model	Maximum run time (d-h:m:s) 3-00:00:00
3D auto-refine	GPU count (0-4) 0 7
3D multi-body	Standard submission script: /net/prog/slurm/relion.sh ? Browse
CIF remement Bayesian polishing	Minimum dedicated cores per node: 1
Mask creation Join star files	Additional arguments:
Particle subtraction	



1. On I/O tab:

One of the 2 unfiltered half-maps: Refine3D/job###/run\_data.star Calibrated pixel size: (from Extract/Polishing job)

2. On ResMap tab: Use ResMap? No

3. On Relion tab: Use Relion? Yes User-provided B-factor: (from Post-process job) MTF of the Detector (STAR file) /net/common/MTF/mtf k2 300kV.star

3. On Running tab: Number of MPI procs: 4 (8 for real data) Number of Threads: 2 Submit to queue? Yes

Click Run! ONCE. This should take ~45 mins.

- 4. In UCSF Chimera:
  - Open relion\_locres\_filtered.mrc and relion\_locres.mrc files from your LocalRes/job### directory
  - 5. Hide the relion\_locres.mrc map
  - On the Volume Viewer, Click on Tools > Surface Color.
  - 7. Color surface relion\_locres\_filtered.mrc by volume data value
  - 8. Set volume file relion\_locres.mrc
  - 9. Click Options
  - 10. Click Reverse
  - 11. Click Color
  - 12. (Optional) Change the histogram contour and click Set (full range of surface values)



#### Cleanup

Relion projects can take up lots of space. In most cases, most of this space can be reclaimed once the processing is done. Relion includes a Gentle clean option to remove the intermediate steps from jobs (we have already used this for removing FOM maps in picking.) This should be run once you're happy that you won't need, for example, iteration 17 from a 25 iteration 2D classification. You can also gently clean an entire project once you've decided that a project has finished.

Relion also includes a Harsh Clean option, which will remove nearly everything apart from the star files that record what you did. You should use this once you've decided a project isn't worth doing anything with, but you still need to keep it, i.e. to match up to your lab notebook, etc,...

Retuber     Retuber       File     Jobs       Schedules     V/O       Reserad pipeline     Alt+R       Catit project note     Alt+R       Print all notes     Alt+D       Overwrite continue     Alt+D       Overwrite continue     Alt+T       About     Alt+Q	<ol> <li>Empty Trash</li> <li>On Relion menu bar, click on File</li> <li>Click on Empty trash</li> </ol>
RELION-3.1-beta-commit-7f8d1d: /home/tjr22/relion31-tutorial	2. Gently clean all jobs
File Jobs Schedules I/O ResMap Relion Running	1 Ón Relion menu har click on Toba
Save job settings Alt+S	
Load job settings Alt+L	2. Click on Gently clean all jobs
Order alphabetically Alt+A	
Under chronologically Alt+L	
Run scheduled inhs	
Stop running scheduled jobs	
Export scheduled job(s)	
Import scheduled job(s)	
Gently clean all jobs Alt+G	
Harshly clean all jobs Alt+H	



#### **Appendix: External jobs**

Relion has the ability to run external jobs and include the results into it's pipeline. So that software written by various groups can be used, we have written a small Python package called relion-external.

The list of currently available packages is:

- Topaz
  - Picking using pre-trained or your own model
  - Training
  - Denoising
- crYOLO
  - Picking using pre-trained or your own model
  - Training
  - Denoising
- Janni
  - Denoising



#### **External: Topaz picking**

Although it is essential that you examine your data (which is why we recommend at least a quick manual-picking strategy,) Al based autopickers perform very well on most datasets. We have integrated Topaz (<u>https://github.com/tbepler/topaz</u>) into Relion via the External job type.

Not all of the features of Topaz are currently available. If you require a feature we haven't yet implements, please let us know. To see what features are available, go to the command line and run:

/net/prog/relion-external/topaz-pick.sh -h

Note that the Relion external job will add "--" to any option label. For example, if you want to specify the model path, use "model" as the label in the Params tab, not "--model". Here we demonstrate how to setup a default run.



Click on the finished job, and display the results in the same was as for any other picking job. You can adjust the default threshold between <u>-6</u> and <u>6</u> (the default is 0). Higher values mean fewer picks. If you add a <u>threshold</u> in the <u>Params</u> tab and continue the job, relion-external will use the already-generated Topaz scores rather than re-running Topaz from scratch.

Note that Topaz generates scores across batches of micrographs - so if you have very poor picking it's possible your problem is not in the micrograph you're looking at but another in the batch. For example: users sometimes accidentally include the gain-reference images output by EPU in their projects. If one of these gain-reference images is included in a batch by Topaz, very poor picking often results.



#### Appendix: 3D classification with multiple initial maps

If you ask Relion to run 3D classification with 4 classes, by default it will make four copies of your initial map, randomly assign each particle to one of the classes, and use the resulting reconstructions as starting points for further classification. The small variations this procedure produce are often more than enough for subsequent iterations to amplify into the classes we want.

If, however, you're concerned that a starting map is too biased towards one particular structure, you can run a 3DAutoRefine job on the whole data set to get a consensus map. That map can then be used as a starting template for further classification.

In rare cases, though, we may want each class to be 'seeded' with very different initial maps (see, for example, Morris et. al, <u>NSMB 2019 26(10), 890-8</u>). We can do this using an appropriately constructed .star file.

For example, if we had four maps in a directory called e.g. relion\_project\_dir/my\_maps/map1.mrc (where relion\_project\_dir/ is our relion project directory,) we could create a star file relion\_project\_dir/my\_maps/initial\_models\_for\_jobXXX.star:

data\_ loop\_ \_rlnReferenceImage my\_maps/map1.mrc my\_maps/map2.mrc my\_maps/map3.mrc my\_maps/map4.mrc

And then put this star file in the Reference Map: field of the IO tab of a 3D Classification job. Note that if you provide 4 maps and ask for more or fewer than 4 classes, Relion will use 4 classes, but subsequent metadata handling will get confused – so make sure that if you specify the same number of classes on the Optimization tab as entries in your .star file.